

F. Heigl,<sup>1,2</sup> M. Murphy,<sup>2</sup> J. Zhou,<sup>2</sup> T.K. Sham,<sup>2</sup> R. Gordon,<sup>3</sup> D. Brewe,<sup>3</sup> S. Heald,<sup>3</sup> H. Shang,<sup>4</sup> and G. Cao<sup>4</sup>


<sup>1</sup> Canadian Synchrotron Radiation Facility, SRC of the University of Wisconsin-Madison

<sup>2</sup> Department of Chemistry, University of Western Ontario, London ON, <sup>3</sup> PNC-CAT, APS Argonne,

<sup>4</sup> Materials Science and Engineering, University of Washington, Seattle




## What is this all about XEOL and timing?



**Timing.** We use the pulsed synchrotron light for a 'pump-probe' like timing experiment. The duration of the pulses are in the order of 100 picoseconds, while the time gap between consecutive bunches of 153 ns gives a relatively wide time window for 'probing' optical photons. The optical decay transitions of many materials have lifetimes which are within this particular time range.

**XEOL.** X-ray excited optical luminescence (XEOL) monitors the decay of an electron-hole pair through optical transitions (infrared to ultraviolet range) after absorption of x-rays. Being sensitive to the size of nanostructures, for example, XEOL can be used to study the phenomenological correlation of optical with structural properties of nanostructures.




## Scintillator Research


Potential scintillator materials emit in the visible light range and have relatively short lifetimes (in the nano second range).

CdWO<sub>4</sub> (solgel) on glass exhibits fast and slow decay components:

- 470nm mainly slow (~56ns)
- 390nm dominates fast (~3ns)





## Nano Materials



**TR-XEOL.** Time resolved XEOL of ZnS:Tb nanowires. Using different time windows shows the time behavior of the fast band gap (mostly within first 20 ns) transition, compared to slow transitions due to defect states.

**Decay.** Lifetime of the XEOL intensity at selected wavelengths (389nm, 497nm, and 532nm).


**Luminescence yield.** XEOL at different excitation energies at the Zn K-edge.




## Metal-Organics

(Ru(bipy)<sub>3</sub>)<sup>2+</sup> is a prototype for transition metal-based photosensitization, charge separation, and photoinduced electron-transfer chemistry. The optical transition is assigned to a metal-to-ligand charge transfer.

### TR-XEOL



### TR-XEOL



Time resolved XEOL reveals: CdWO<sub>4</sub> fast and slow emission lines in the visible light range, after excitation at 10213eV photon energy.

## Acknowledgement

We greatly acknowledge Richard Rosenberg for help and technical advice, as well the following funding agencies and institutions:

